

# How to use your Anapen®

For emergency treatment of anaphylaxis

## Anapen®



500 micrograms adrenaline (epinephrine) in 0.3ml solution for injection in a pre-filled syringe.  
For adults weighing over 65kg or at high risk of anaphylaxis.

## Anapen®



300 micrograms adrenaline (epinephrine) in 0.3 ml solution for injection in a pre-filled syringe.  
For adults weighing under 65kg and children, depending on the body weight of the child and the discretion of the prescriber.

## Anapen® junior



150 micrograms adrenaline (epinephrine) in 0.3 ml solution for injection in a pre-filled syringe.  
For children, depending on the body weight of the child and the discretion of the prescriber

For further information please refer to the patient information leaflet.

To order Anapen please contact:

### BOC Healthcare

Customer Service Centre,  
Priestley Road, Worsley,  
Manchester, M28 2UT  
Phone: 0161 930 6010  
Email: [bochealthcare-uk@boc.com](mailto:bochealthcare-uk@boc.com)  
Website: [www.boclifeline.co.uk](http://www.boclifeline.co.uk)

For clinical support please contact:

Phone: +44 (0) 1748 828785  
Email: [info@lincolnmedical.co.uk](mailto:info@lincolnmedical.co.uk)  
Website: [www.anapen.co.uk](http://www.anapen.co.uk)

Please visit **[www.anapen.co.uk](http://www.anapen.co.uk)** to see a video and an animation on "how to use Anapen" and for full Summary of Product Characteristics documents and Reminder Registration facilities.

# How to use your Anapen<sup>®</sup>



1. Remove the black needle cap



2. Remove the safety cap from the red firing button



3. Hold Anapen<sup>®</sup> against the outer thigh and press the red firing button



4. Hold Anapen<sup>®</sup> in position for 10 seconds to allow the full dose of adrenaline to be injected. Gently massage the injection site