

ENTONOX[®] demand valve

Maximum performance and protection

The ENTONOX demand valve delivers the best available gas flow performance with the most comfortable handset on the market whilst having a long warranty and servicing interval designed specifically for the demanding acute healthcare environment.

Features and benefits

- Very low resistance for reduced patient effort.
- Patented 'exhalation valve' prevents expired breath entering the handset.
- Single patient use exhalation valve with filter prevents cross contamination.
- Rugged design with 4 year warranty.
- Compact, lightweight design for improved patient comfort.
- Only approved polymers in contact with the high pressure gas (as per ISO 15001).
- Wipe clean soft cover and strap.
- 4 year service interval.

Assured quality/standards

All Carnét demand valves have been tested to the exacting technical standards, as specified by the latest ISO standard, allowing them to be CE Marked to the Medical Device Directive, giving you the confidence that comes with a safe and reliable medical device. Carnét demand valves meet or surpass the requirements of:

- BS 4272-2
- BS 5682
- BS EN 1281-1
- BS EN 739
- BS EN ISO 15001




Specifications

Inspiratory resistance	< 1.5 kPa at 200 l/min < 0.25 kPa at 10 l/min
Peak flow	> 200 l/min
Response time	< 100 ms to 200 l/min
Environment	Storage: -20°C to + 60°C Use: 0°C to + 40°C
CE marking	Active Medical Device - Class IIa
Warranty period	4 years
Service interval	4 years
Product order codes	Specified overleaf

Demand valve

Order code	Description
828-0013	Demand valve – ENTONOX – 2m hose
828-0001	Demand valve – ENTONOX – 3m hose
828-0014	Demand valve – ENTONOX – 4m hose
828-0015	Demand valve – ENTONOX – 5m hose
828-0016	Demand valve – ENTONOX – 6m hose
828-0002	Exhalation valve – single patient use – 100 pack
828-0019	Exhalation valve – AGSS – single patient use – 50 pack
212-0067	Demand valve – Lanyard



To find out more or place your order today contact:

08456 447 764

or for 24 hour support visit:

www.bochealthcare.co.uk

BOC Healthcare, Great Britain & N. Ireland
Customer Service Centre, Priestley Road,
Worsley, Manchester M28 2UT, United Kingdom
Tel 08456 447 764, Fax 0845 087 0062
bochealthcare-uk@boc.com, www.bochealthcare.co.uk

BOC Healthcare, Republic of Ireland
PO Box 201, Bluebell, Dublin 12, Republic of Ireland
Tel 01 409 1827, Fax 01 409 1854
healthcareinfo.ie@boc.com, www.BOCOnline.ie