

12 cylinder trolley

For easy multiple cylinder transportation

The 12 cylinder trolley from BOC Healthcare can accommodate up to 12 medical gas cylinders for safe and efficient transportation.

Ergonomically designed with a modern healthcare environment in mind, this new design from BOC Healthcare will allow the safe movement of multiple cylinders. The design details will benefit portering and clinical staff.

Features and benefits

- Twin handled design makes it easy to manoeuvre the trolley
- End plates contain signage for HTM-02 compliance
- Rubber base mat aids stability and reduces noise by stopping cylinder movement
- Wheel bumpers help reduce upper trolley impacts
- Tray heights appropriate for secure support of various 100mm diameter cylinders
- Retaining trays contoured to aid ease of loading and unloading and minimise damage to cylinders
- 4 x 125 mm diameter wheels
- 2 x lockable wheel brakes
- Designed to accommodate 100mm diameter cylinders, e.g. CD, D, ZD, ED, and E sizes

Important notes

- Ensure cylinders are loaded evenly to help spread the load and to increase stability and manoeuvrability
- Always assess whether it is safe to move a trolley across ground that is uneven
- Do not attempt to load more than 12 cylinders

Product code: 19321752



To find out more contact
08456 447 764 24 hours a day, or visit
www.bochealthcare.co.uk

Specifications

Handle height	1063 mm	Width	450mm (520mm with wheel bumpers)
Upper tray height	550 mm	Length	895 mm
Wheels	125 mm diameter rubber	Finish	Powder coated
Bumpers	Neoprene	Cleaning	Soap and water with soft cloth or brush
Weight	30 kg (unladen)	Disinfecting	Virkon solution
Maximum load	12 cylinders loaded into tray holes only	Labelling	2 cylinder storage labels
Frame	Steel		

Product code: 19321752



BOC Healthcare, Great Britain & N. Ireland
Customer Service Centre, Priestley Road,
Worsley, Manchester M28 2UT, United Kingdom
Tel 08456 447 764, Fax 0845 087 0062
bochealthcare-uk@boc.com, www.bochealthcare.co.uk

BOC Healthcare, Republic of Ireland
PO Box 201, Bluebell, Dublin 12, Republic of Ireland
Tel 01 409 1827, Fax 01 409 1854
healthcareinfo.ie@boc.com, www.BOConline.ie